## Structure of the product and materials

Polyphenylsulfone body – PPSU black  
Invasive elements in steel pursuant to ASTM F138 – ISO 5832-2

## Method of installation

### Proximal interphalangeal joint

Placement of a first d: 1.5 mm Schanz pin using a drill at low speed.

Percutaneous insertion of the first pin is made in the neck of the proximal phalanx through both cortices. The percutaneous insertion of the second Schanz Pin is made on the base of the middle phalanx through both cortices.

Percutaneous insertion of one or two supplementary pins is made in the diaphyses of the proximal and middle phalanges through both cortices.

**The Digifix fixator is installed and operational.**

### Metacarpophalangeal joint

Placement of a first d: 1.5 mm Schanz pin using a drill at low speed.

Percutaneous insertion of the first pin is made in the neck of the metacarpal phalanx through both cortices. The percutaneous insertion of the second Schanz Pin is made in the base of the proximal phalanx through both cortices.

Percutaneous insertion of one or two supplementary pins is made in the diaphyses of the metacarpal and middle phalanges through both cortices.

**The Digifix fixator is installed and operational.**

## Advantages and benefits

- > Pivot lock- Early mobilization of the joint
- > Setting the relaxation force
- > Composite materials – light for the patient
- > Radio transparent, allows better monitoring of the trauma
- > Two possible positions for each pin
- > Exercise protocol available for each patient

## Indications / use

- > IFP Articular fractures
- > MF Articular fractures
- > Open fractures of phalanges and/or associated with loss of soft tissue
- > Unstable metaphyseal fractures

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