

AMPLIFIX® 1 and 3 AMPLIFIX® 1 and 3 with Gentamicin

Surgical cement



References



AMPLIFIX 1 and 3

Amplifix 1 - Surgical cement - Standard Viscosity	1-0400101
Amplifix 3 - Surgical cement - Low Viscosity	1-0400103

Indication: AMPLIFIX 1&3 bone cement is indicated for fixation of prostheses to living bone in orthopaedic musculoskeletal surgical procedures for rheumatoid arthritis, osteoarthritis, traumatic arthritis, osteoporosis, avascular necrosis, collagen disease, severe joint destruction secondary to trauma or other conditions, and revision of previous arthroplasty. MD class : IIb

AMPLIFIX 1 and 3 - Gentamicin

Amplifix 1 G - Surgical cement with gentamicin - Standard Viscosity	1-0400201
Amplifix 3 G - Surgical cement with gentamicin - Low Viscosity	1-0400203

Indication: AMPLIFIX bone cements with gentamicin are indicated for use as bone cement in arthroplasty procedures of the hip, knee and to fix plastic or metal prosthetic parts to bone. These bone cements are indicated for prostheses fixation in primary or secondary surgeries as well as to affix a new prosthesis in the second phase of a two-stage revision after the initial infection has been cleared. MD class: III

Packaging

MONOMER : sterile liquid ampoule in a blister

POLYMER : sterile bag powder (40 grams for the version without gentamicin and 41,6 grams for the version with gentamicin)

Double sterile packaging

STORAGE : maximum temperature of 25 °C

Legal Manufacturer: TEKNIMED SAS, 8, rue du Corps Franc-Pommiès, 65500 VIC en BIGORRE, France

Distributed by AMPLITUDE

Service Clients-France :

Porte du Grand Lyon,
01700 Neyron, France
Tel. : +33 (0)4 37 85 19 19
Fax : +33 (0)4 37 85 19 18
E-mail : amplitude@amplitude-ortho.com

Customer Service-Export :

11, cours Jacques Offenbach. ZA Mozart 2,
26000 Valence, France
Tel. : +33 (0)4 75 41 87 41
Fax : +33 (0)4 75 41 87 42
E-mail : amplitude@amplitude-ortho.com



AMPLITUDE[®]

Internet : www.amplitude-ortho.com